

Product Datasheet

Mouse Fibroblast Growth Factor Receptor 3 (FGFR3) ELISA Kit (orb1208524)

Description	Mouse Fibroblast Growth Factor Receptor 3 (FGFR3) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q61851
Sensitivity	0.057ng/mL
Entrez	14184
Expiration Date	Please enquire.