

Product Datasheet

Mouse Fibroblast Growth Factor Receptor 1 (FGFR1) ELISA Kit (orb1208525)

Description	Mouse Fibroblast Growth Factor Receptor 1 (FGFR1) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	P16092
Sensitivity	0.061ng/mL
Entrez	14182
Expiration Date	Please enquire.