

Product Datasheet

Rabbit Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit (orb1208532)

Description	Rabbit Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Rabbit
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48799
Sensitivity	6.4pg/mL
Entrez	100009068
Expiration Date	Please enquire.