

Product Datasheet

Mouse Epidermal Growth Factor Receptor 2 (EGFR2) ELISA Kit (orb1208540)

Description	Mouse Epidermal Growth Factor Receptor 2 (EGFR2) ELISA Kit
Reactivity	Mouse
Range	0.625-40ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P70424
Sensitivity	0.207ng/mL
Entrez	13866
Expiration Date	Please enquire.