

Product Datasheet

Rat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit (orb1208547)

Description	Rat Bone Morphogenetic Protein 4 (BMP4) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06826
Sensitivity	6.1pg/mL
Entrez	25296
Expiration Date	Please enquire.