

Product Datasheet

Rabbit Brain Derived Neurotrophic Factor (BDNF) ELISA Kit (orb1208552)

Description	Rabbit Brain Derived Neurotrophic Factor (BDNF) ELISA Kit
Reactivity	Rabbit
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	G1SRK9
Sensitivity	0.055ng/mL
Entrez	100356949
Expiration Date	Please enquire.