

Product Datasheet

Rat Adiponectin Receptor 2 (ADIPOR2) ELISA Kit (orb1208564)

Description	Rat Adiponectin Receptor 2 (ADIPOR2) ELISA Kit
Reactivity	Rat
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates or other biological fluids.
Uniprot ID	Q3YAF7
Sensitivity	0.125ng/mL
Entrez	312670
Expiration Date	Please enquire.