

Product Datasheet

Mouse Fatty Acid Binding Protein 2, Intestinal (FABP2) ELISA Kit (orb1208577)

Description	Mouse Fatty Acid Binding Protein 2, Intestinal (FABP2) ELISA Kit
Reactivity	Mouse
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P55050
Sensitivity	0.063ng/mL
Entrez	14079
Expiration Date	Please enquire.