

Product Datasheet

Bovine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit (orb1208589)

Description	Bovine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Bovine
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06599
Sensitivity	32pg/mL
Entrez	280943
Expiration Date	Please enquire.