

Product Datasheet

Human Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit (orb1208596)

Description	Human Tumor Necrosis Factor Ligand Superfamily, Member 7 (TNFSF7) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P32970
Sensitivity	0.056ng/mL
Entrez	970
Expiration Date	Please enquire.