

Product Datasheet

Human Milk Fat Globule EGF Factor 8 (MFGE8) ELISA Kit (orb1208601)

Description	Human Milk Fat Globule EGF Factor 8 (MFGE8) ELISA Kit
Reactivity	Human
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	Q08431
Sensitivity	12.3pg/mL
Entrez	4240
Expiration Date	Please enquire.