

Product Datasheet

Human Latent Transforming Growth Factor Beta Binding Protein 2 (LTBP2) ELISA Kit (orb1208603)

Description	Human Latent Transforming Growth Factor Beta Binding Protein 2 (LTBP2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q14767
Sensitivity	0.062ng/mL
Entrez	4053
Expiration Date	Please enquire.