

Product Datasheet

Human Lymphocyte Function Associated Antigen 1 Alpha (LFA1a) ELISA Kit (orb1208605)

Description	Human Lymphocyte Function Associated Antigen 1 Alpha (LFA1a) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P20701
Sensitivity	0.055ng/mL
Entrez	3683
Expiration Date	Please enquire.