

## Product Datasheet

# Human Interferon Regulatory Factor 5 (IRF5) ELISA Kit (orb1208607)

|                          |                                                               |
|--------------------------|---------------------------------------------------------------|
| <b>Description</b>       | Human Interferon Regulatory Factor 5 (IRF5) ELISA Kit         |
| <b>Reactivity</b>        | Human                                                         |
| <b>Range</b>             | 0.156-10ng/mL                                                 |
| <b>Note</b>              | For research use only                                         |
| <b>Application notes</b> | sandwich                                                      |
| <b>Sample Types</b>      | serum, plasma, tissue homogenates or other biological fluids. |
| <b>Uniprot ID</b>        | <b>Q13568</b>                                                 |
| <b>Sensitivity</b>       | 0.051ng/mL                                                    |
| <b>Entrez</b>            | <b>3663</b>                                                   |
| <b>Expiration Date</b>   | Please enquire.                                               |