

Product Datasheet

Human Interleukin 1 Receptor Associated Kinase 4 (IRAK4) ELISA Kit (orb1208610)

Description	Human Interleukin 1 Receptor Associated Kinase 4 (IRAK4) ELISA Kit
Reactivity	Human
Range	1.56-100ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q9NWZ3
Sensitivity	0.53ng/mL
Entrez	51135
Expiration Date	Please enquire.