

Product Datasheet

Human Interleukin 12 Receptor Beta 2 (IL12Rb2) ELISA Kit (orb1208634)

Description	Human Interleukin 12 Receptor Beta 2 (IL12Rb2) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	Q99665
Sensitivity	0.055ng/mL
Entrez	3595
Expiration Date	Please enquire.