

Product Datasheet

Human Interleukin 12 Receptor Beta 1 (IL12Rb1) ELISA Kit (orb1208635)

Description	Human Interleukin 12 Receptor Beta 1 (IL12Rb1) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P42701
Sensitivity	0.071ng/mL
Entrez	3594
Expiration Date	Please enquire.