

Product Datasheet

Human Fibroblast Growth Factor 10 (FGF10) ELISA Kit (orb1208646)

Description	Human Fibroblast Growth Factor 10 (FGF10) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O15520
Sensitivity	5.7pg/mL
Entrez	2255
Expiration Date	Please enquire.