

Product Datasheet

Human Cluster Of Differentiation 8a (CD8a) ELISA Kit (orb1208686)

Description	Human Cluster Of Differentiation 8a (CD8a) ELISA Kit
Reactivity	Human
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01732
Sensitivity	1.3ng/mL
Entrez	925
Expiration Date	Please enquire.