

Product Datasheet

Mouse Transforming Growth Factor Alpha (TGFA) ELISA Kit (orb1208711)

Description	Mouse Transforming Growth Factor Alpha (TGFA) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P48030
Sensitivity	5.5pg/mL
Entrez	21802
Expiration Date	Please enquire.