

Product Datasheet

Mouse Platelet Derived Growth Factor Subunit A (PDGFA) ELISA Kit (orb1208712)

Description	Mouse Platelet Derived Growth Factor Subunit A (PDGFA) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	P20033
Sensitivity	12.9pg/mL
Entrez	18590
Expiration Date	Please enquire.