

Product Datasheet

Mouse Interleukin 1 Receptor Type II (IL1R2) ELISA Kit (orb1208719)

Description	Mouse Interleukin 1 Receptor Type II (IL1R2) ELISA Kit
Reactivity	Mouse
Range	31.2-2000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P27931
Sensitivity	11.38pg/mL
Entrez	16178
Expiration Date	Please enquire.