

Product Datasheet

Mouse Interleukin 1 Receptor Type I (IL1R1) ELISA Kit (orb1208720)

Description	Mouse Interleukin 1 Receptor Type I (IL1R1) ELISA Kit
Reactivity	Mouse
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates or other biological fluids.
Uniprot ID	P13504
Sensitivity	29pg/mL
Entrez	16177
Expiration Date	Please enquire.