

Product Datasheet

Mouse Insulin Like Growth Factor Binding Protein 4 (IGFBP4) ELISA Kit (orb1208732)

Description	Mouse Insulin Like Growth Factor Binding Protein 4 (IGFBP4) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P47879
Sensitivity	6.4pg/mL
Entrez	16010
Expiration Date	Please enquire.