

Product Datasheet

Mouse Fibroblast Growth Factor 9 (FGF9) ELISA Kit (orb1208740)

Description	Mouse Fibroblast Growth Factor 9 (FGF9) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P54130
Sensitivity	6.6pg/mL
Entrez	14180
Expiration Date	Please enquire.