

Product Datasheet

Mouse Adiponectin Receptor 1 (ADIPOR1) ELISA Kit (orb1208760)

Description	Mouse Adiponectin Receptor 1 (ADIPOR1) ELISA Kit
Reactivity	Mouse
Range	3.12-200ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma or other biological fluids.
Uniprot ID	Q91VH1
Sensitivity	1.17ng/mL
Entrez	72674
Expiration Date	Please enquire.