

Product Datasheet

Human Lymphocyte Function Associated Antigen 3 (LFA3) ELISA Kit (orb1208781)

Description	Human Lymphocyte Function Associated Antigen 3 (LFA3) ELISA Kit
Reactivity	Human
Range	12.5-800pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, cell lysates or other biological fluids.
Uniprot ID	P19256
Sensitivity	5.8pg/mL
Entrez	965
Expiration Date	Please enquire.