

Product Datasheet

Human Interleukin 8 Receptor Alpha (IL8Ra) ELISA Kit (orb1208784)

Description	Human Interleukin 8 Receptor Alpha (IL8Ra) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P25024
Sensitivity	0.072ng/mL
Entrez	3577
Expiration Date	Please enquire.