

Product Datasheet

Human Interferon Alpha 4 (IFN α 4) ELISA Kit (orb1208802)

Description	Human Interferon Alpha 4 (IFN α 4) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P05014
Sensitivity	5.4pg/mL
Entrez	3441
Expiration Date	Please enquire.