

Product Datasheet

Human Interleukin 1 Delta (FIL1d) ELISA Kit (orb1208805)

Description	Human Interleukin 1 Delta (FIL1d) ELISA Kit
Reactivity	Human
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9UBH0
Sensitivity	5.7pg/mL
Entrez	26525
Expiration Date	Please enquire.