

Product Datasheet

Human Bone Morphogenetic Protein Receptor 1A (BMPR1A) ELISA Kit (orb1208813)

Description	Human Bone Morphogenetic Protein Receptor 1A (BMPR1A) ELISA Kit
Reactivity	Human
Range	0.312-20ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	P36894
Sensitivity	0.11ng/mL
Entrez	657
Expiration Date	Please enquire.