

Product Datasheet

Human Adiponectin Receptor 2 (ADIPOR2) ELISA Kit (orb1208822)

Description	Human Adiponectin Receptor 2 (ADIPOR2) ELISA Kit
Reactivity	Human
Range	0.781-50ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	tissue homogenates, cell lysates or other biological fluids.
Uniprot ID	Q86V24
Sensitivity	0.24ng/mL
Entrez	79602
Expiration Date	Please enquire.