

Product Datasheet

Rat Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit (orb1208834)

Description	Rat Fibroblast Growth Factor 2, Basic (FGF2) ELISA Kit
Reactivity	Rat
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P13109
Sensitivity	4.37pg/mL
Entrez	54250
Expiration Date	Please enquire.