

Product Datasheet

Mouse Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) ELISA Kit (orb1208861)

Description	Mouse Heparin Binding Epidermal Growth Factor Like Growth Factor (HBEGF) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q06186
Sensitivity	6.0pg/mL
Entrez	15200
Expiration Date	Please enquire.