

Product Datasheet

Canine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit (orb1208878)

Description	Canine Tumor Necrosis Factor Alpha (TNFa) ELISA Kit
Reactivity	Canine
Range	3.12-200pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P51742
Sensitivity	1.19pg/mL
Entrez	403922
Expiration Date	Please enquire.