

Product Datasheet

Human Kallikrein 10 (KLK10) ELISA Kit (orb1208882)

Description	Human Kallikrein 10 (KLK10) ELISA Kit
Reactivity	Human
Range	0.156-10ng/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	O43240
Sensitivity	0.048ng/mL
Entrez	5655
Expiration Date	Please enquire.