

Product Datasheet

Guinea pig Interleukin 8 (IL8) ELISA Kit (orb1208883)

Description	Guinea pig Interleukin 8 (IL8) ELISA Kit
Reactivity	Guinea pig
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P49113
Sensitivity	6.8pg/mL
Entrez	100379599
Expiration Date	Please enquire.