

Product Datasheet

Rabbit Interleukin 1 Beta (IL1b) ELISA Kit (orb1208965)

Description	Rabbit Interleukin 1 Beta (IL1b) ELISA Kit
Reactivity	Rabbit
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P14628
Sensitivity	33pg/mL
Entrez	100008990
Expiration Date	Please enquire.