

Product Datasheet

Mouse Interferon Alpha (IFNa) ELISA Kit (orb1208987)

Description	Mouse Interferon Alpha (IFNa) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	P01572
Sensitivity	6.3pg/mL
Entrez	111654
Expiration Date	Please enquire.