

Product Datasheet

Human Interleukin 1 Zeta (IL1z) ELISA Kit (orb1208996)

Description	Human Interleukin 1 Zeta (IL1z) ELISA Kit
Reactivity	Human
Range	78.1-5000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9NZH6
Sensitivity	24pg/mL
Entrez	27178
Expiration Date	Please enquire.