

Product Datasheet

Canine Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit (orb1209024)

Description	Canine Vascular Endothelial Growth Factor A (VEGFA) ELISA Kit
Reactivity	Canine
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9MYV3
Sensitivity	5.4pg/mL
Entrez	403802
Expiration Date	Please enquire.