

Product Datasheet

Mouse Interleukin 23 (IL23) ELISA Kit (orb1209044)

Description	Mouse Interleukin 23 (IL23) ELISA Kit
Reactivity	Mouse
Range	15.6-1000pg/mL
Note	For research use only
Application notes	sandwich
Sample Types	serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Uniprot ID	Q9EQ14
Sensitivity	5.1pg/mL
Entrez	83430
Expiration Date	Please enquire.