

Product Datasheet

Human 5-hydroxytryptamine ELISA Kit (orb1209120)

Description	Human 5-hydroxytryptamine, 5HT ELISA Kit
Reactivity	Human
Range	7 - 1500ng/ml
Tested Applications	ELISA
Target	5HT
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	1600ng/mL
Sensitivity	4.33ng/ml
Expiration Date	Please enquire.