

Product Datasheet

Fish Superoxide Dismutase ELISA Kit (orb1209128)

Description	Fish Superoxide Dismutase, SOD ELISA Kit
Reactivity	Fish
Range	0.1 - 40ng/ml
Tested Applications	ELISA
Target	SOD
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	48ng/mL
Sensitivity	0.045ng/ml
Expiration Date	Please enquire.