

Product Datasheet

Fish Glutathione Peroxidase ELISA Kit (orb1209130)

Description	Fish Glutathione Peroxidase, GSH-PX ELISA Kit
Reactivity	Fish
Range	50 - 10000uU/ml
Tested Applications	ELISA
Target	GSH-Px
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	12800μU/mL
Sensitivity	23.95μU/ml
Expiration Date	Please enquire.