

Product Datasheet

Bovine Adiponectin ELISA Kit (orb1209134)

Description	Bovine Adiponectin, ADP ELISA Kit
Reactivity	Bovine
Range	0.3 - 90ug/ml
Tested Applications	ELISA
Target	ADP
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	96ug/ml
Sensitivity	0.17µg/ml
Expiration Date	Please enquire.