

Product Datasheet

Sheep Progesterone ELISA Kit (orb1209142)

Description	Sheep Progesterone, PROG ELISA Kit
Reactivity	Sheep
Range	0.2 - 60ng/ml
Tested Applications	ELISA
Target	PROG
Storage	2-8°C
Note	For research use only
Application notes	Competitive
Sample Types	Serum, plasma, cell culture supernatant
Solubility (25°C)	64ng/ml
Sensitivity	0.097ng/ml
Expiration Date	Please enquire.