

## Product Datasheet

### Bovine Growth Hormone ELISA Kit (orb1209147)

|                            |                                         |
|----------------------------|-----------------------------------------|
| <b>Description</b>         | Bovine Growth Hormone,<br>GH ELISA Kit  |
| <b>Reactivity</b>          | Bovine                                  |
| <b>Range</b>               | 0.1 - 40ng/ml                           |
| <b>Tested Applications</b> | ELISA                                   |
| <b>Target</b>              | GH                                      |
| <b>Storage</b>             | 2-8°C                                   |
| <b>Note</b>                | For research use only                   |
| <b>Application notes</b>   | Competitive                             |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernatant |
| <b>Solubility (25°C)</b>   | 48ng/ml                                 |
| <b>Sensitivity</b>         | 0.046ng/ml                              |
| <b>Expiration Date</b>     | Please enquire.                         |