

Product Datasheet

Human Mitogen-activated Protein Kinase 4 ELISA Kit (orb1209149)

Description	Human Mitogen-activated Protein Kinase 4, MAP4K4 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	MAP4K4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	142,101
Uniprot ID	O95819
Solubility (25°C)	3200ng/L
Sensitivity	7.83ng/L
Entrez	9448
Expiration Date	Please enquire.