

Product Datasheet

Human Immunoglobulin Lambda-like Polypeptide 1 ELISA Kit (orb1209154)

Description	Human Immunoglobulin Lambda-like Polypeptide 1, IGLL1 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	IGLL1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	22,963
Uniprot ID	P15814
Solubility (25°C)	2400ng/L
Sensitivity	5.26ng/L
Entrez	3543
Expiration Date	Please enquire.